

BIOPSY FORCEPS


disposable



ENDO-FLEX[®]


disposable Biopsy Forceps \varnothing 2.3 mm, for working channels $\geq \varnothing$ 2.8 mm

stainless steel spiral, plastic handle:



length	package	oval fenestrated cups	
		without spike	with spike
160 cm	10 pieces	NE 02 22-B	NE 02 22D-B
180 cm	10 pieces	NE 02 22-C	NE 02 22D-C
230 cm	10 pieces	NE 02 22-G	NE 02 22D-G


coated stainless steel spiral, plastic handle:



length	package	oval fenestrated cups	
		without spike	with spike
160 cm	10 pieces	NEX 02 22-B	NEX 02 22D-B
230 cm	10 pieces	NEX 02 22-G	NEX 02 22D-G

disposable Biopsy Forceps \varnothing 1.8 mm, for Bronchoscopy, for working channels $\geq \varnothing$ 2.0 mm

stainless steel spiral, plastic handle:



length	package	oval fenestrated cups without spike
120 cm	10 pieces	NE 16 18-A

Legend:



This instrument is designed and warranted for single use only.



This instrument is sterile packed separately.

ENDO-FLEX GmbH
Alte Hünxer Str. 115
46562 Voerde / Germany
Phone: +49 281 - 9 44 00 - 0
Fax: +49 281 - 9 44 00 - 11
e-mail: export@endo-flex.de
www.endo-flex.de