


POLYPECTOMY SNARES

disposable



 **ENDO-FLEX**[®]

rotatable disposable Polypectomy Snare \varnothing 2.3 mm, plastic handle, for working channels $\geq \varnothing$ 2.8 mm



STERILE EO



oval loop

braided wire

length	packing unit	15 mm	25 mm	35 mm
180 cm	5 pieces	NOE 34 2214-C	NOE 34 2216-C	NOE 34 2217-C
230 cm	5 pieces	NOE 34 2214-G	NOE 34 2216-G	NOE 34 2217-G



STERILE EO



asymmetrical loop

braided wire

length	packing unit	15 mm	25 mm
230 cm	5 pieces	NOE 3222-G	NOE 3322-G



STERILE EO



oval loop

monofilament wire

length	packing unit	15 mm	25 mm
230 cm	5 pieces	NOE 34 2214M-G	NOE 34 2216M-G

recommended equipment

active cords:



HF generator	Order No.
ERBE E- and T-Series	64 0300
ERBE ICC	64 0500
OLYMPUS HF 120/300	64 0500



Legend:



This instrument is designed and warranted for single use only.



This instrument is designed to be used with HF current.



This instrument is rotatable.

STERILE EO

This instrument is sterile packed separately.

ENDO-FLEX GmbH
Alte Hünxer Str. 115
46562 Voerde / Germany
Phone: +49 281 - 9 44 00 - 0
Fax: +49 281 - 9 44 00 - 11
e-mail: export@endo-flex.de
www.endo-flex.de