


Polypectomy Snare \varnothing 2.3 mm, incl. handle,
for working channels $\geq \varnothing$ 2.8 mm

asymmetrical loop



length	40 mm		60 mm	
	standard tip	insulated tip	standard tip	insulated tip
180 cm	OK 65 2240-C	OK 65 22401-C	OK 65 2260-C	OK 65 22601-C
230 cm	OK 65 2240-G	OK 65 22401-G	OK 65 2260-G	OK 65 22601-G

Polypectomy Snare \varnothing 2.3 mm, fixed handle,
for working channels $\geq \varnothing$ 2.8 mm

braided wire
thermo-coated
standard tip
oval loop



braided wire
semi-rigid
standard tip
oval loop



length	35 mm	20 - 35 mm
230 cm	OK 34 2237-G	OK 34 2267-G

recommended equipment

active cords:



HF generators	Order No.
ERBE E- and T-Series	64 0300
ERBE ICC	64 0500
OLYMPUS HF 120/300	64 0500



universal handle:

(for snares, which are supplied without handle)

Order No.
62 0300



Legend:



This instrument is autoclavable by 134°C / 273°F.



This instrument is detachable for cleaning.



This instrument is designed to be used with HF current.




This Instrument possesses a flushing port.

ENDO-FLEX GmbH
Alte Hünxer Str. 115
46562 Voerde / Germany
Phone: +49 281 - 9 44 00 - 0
Fax: +49 281 - 9 44 00 - 11
e-mail: export@endo-flex.de
www.endo-flex.de


© ENDO-FLEX - Technical specifications are subject to change

POLYPECTOMY SNARES
re-usable



Polypectomy Snare \varnothing 2.3 mm, fixed handle, for working channels $\geq \varnothing$ 2.8 mm

asymmetrical loop



length	15 mm		25 mm	
	insulated tip	standard tip	insulated tip	standard tip
180 cm	OK 30 22-C	OK 32 22-C	OK 31 22-C	OK 33 22-C
230 cm	OK 30 22-G	OK 32 22-G	OK 31 22-G	OK 33 22-G

oval loop



length	20 mm		25 mm		35 mm	
	insulated tip	standard tip	insulated tip	standard tip	insulated tip	standard tip
180 cm	OK 34 2205-C	OK 34 2215-C	OK 34 2206-C	OK 34 2216-C	OK 34 2207-C	OK 34 2217-C
230 cm	OK 34 2205-G	OK 34 2215-G	OK 34 2206-G	OK 34 2216-G	OK 34 2207-G	OK 34 2217-G

hexagonal loop



length	20 mm		25 mm		35 mm	
	insulated tip	standard tip	insulated tip	standard tip	insulated tip	standard tip
180 cm	OK 38 2205-C	OK 37 2205-C	OK 38 2206-C	OK 37 2206-C	OK 38 2207-C	OK 37 2207-C
230 cm	OK 38 2205-G	OK 37 2205-G	OK 38 2206-G	OK 37 2206-G	OK 38 2207-G	OK 37 2207-G

oval loop monofilament wire



length	20 mm	25 mm	35 mm
	standard tip	standard tip	standard tip
180 cm	OK 34 2215M-C	OK 34 2216M-C	OK 34 2217M-C
230 cm	OK 34 2215M-G	OK 34 2216M-G	OK 34 2217M-G

Polypectomy Snare \varnothing 2.3 mm, without handle, for working channels $\geq \varnothing$ 2.8 mm

asymmetrical loop



length	15 mm		25 mm		15 mm		25 mm	
	insulated tip				standard tip			
	standard	rotatable	standard	rotatable	standard	rotatable	standard	rotatable
180 cm	30 22-C	D 30 22-C	31 22-C	D 31 22-C	32 22-C	D 32 22-C	33 22-C	D 33 22-C
230 cm	30 22-G	D 30 22-G	31 22-G	D 31 22-G	32 22-G	D 32 22-G	33 22-G	D 33 22-G
280 cm							33 22-M	
300 cm							33 22-O	

insulated tip oval loop



length	20 mm		25 mm		35 mm	
	standard	rotatable	standard	rotatable	standard	rotatable
180 cm	34 2205-C	D 34 2205-C	34 2206-C	D 34 2206-C	34 2207-C	D 34 2207-C
230 cm	34 2205-G	D 34 2205-G	34 2206-G	D 34 2206-G	34 2207-G	D 34 2207-G
280 cm			34 2206-M			
300 cm			34 2206-O			

standard tip oval loop



length	20 mm		25 mm		35 mm	
	standard	rotatable	standard	rotatable	standard	rotatable
180 cm	34 2215-C	D 34 2215-C	34 2216-C	D 34 2216-C	34 2217-C	D 34 2217-C
230 cm	34 2215-G	D 34 2215-G	34 2216-G	D 34 2216-G	34 2217-G	D 34 2217-G
280 cm	34 2215-M		34 2216-M			
300 cm	34 2215-O		34 2216-O			

insulated tip hexagonal loop



length	20 mm		25 mm		35 mm	
	standard	rotatable	standard	rotatable	standard	rotatable
180 cm	38 2205-C	D 38 2205-C	38 2206-C	D 38 2206-C	38 2207-C	D 38 2207-C
230 cm	38 2205-G	D 38 2205-G	38 2206-G	D 38 2206-G	38 2207-G	D 38 2207-G

necessary devices: universal handle 62 0300



Polypectomy Snare \varnothing 2.3 mm, without handle, for working channels $\geq \varnothing$ 2.8 mm

standard tip hexagonal loop



length	20 mm		25 mm		35 mm	
	standard	rotatable	standard	rotatable	standard	rotatable
180 cm	37 2205-C	D 37 2205-C	37 2206-C	D 37 2206-C	37 2207-C	D 37 2207-C
230 cm	37 2205-G	D 37 2205-G	37 2206-G	D 37 2206-G	37 2207-G	D 37 2207-G

standard tip oval loop

monofilament wire



length	20 mm		25 mm		35 mm	
	standard	rotatable	standard	rotatable	standard	rotatable
180 cm	34 2215M-C	D 34 2215M-C	34 2216M-C	D 34 2216M-C	34 2217M-C	D 34 2217M-C
230 cm	34 2215M-G	D 34 2215M-G	34 2216M-G	D 34 2216M-G	34 2217M-G	D 34 2217M-G

Polypectomy Snare \varnothing 1.8 mm, without handle, for working channels $\geq \varnothing$ 2.0 mm

oval loop



length	20 mm		25 mm	20 mm		25 mm
	insulated tip		insulated tip	standard tip		standard tip
180 cm	34 1805-C		34 1806-C	34 1815-C		34 1816-C
230 cm	34 1805-G		34 1806-G	34 1815-G		34 1816-G
280 cm				34 1815-M		
300 cm				34 1815-O		



length	15 mm	25 mm	20 mm	25 mm	20 mm	25 mm
	asymmetrical loop, standard tip		hexagonal loop, insulated tip		hexagonal loop, standard tip	
180 cm	32 18-C	33 18-C	38 1805-C	38 1806-C	37 1805-C	37 1806-C
230 cm	32 18-G	33 18-G	38 1805-G	38 1806-G	37 1805-G	37 1806-G

necessary devices: universal handle 62 0300

