

# SCLEROTHERAPY NEEDLES

disposable



© ENDO-FLEX - Technical specification are subject to change

## Legend:



This instrument is designed and warranted for single use only.



This instrument is sterile packed separately.

ENDO-FLEX GmbH  
Alte Hünxer Str. 115  
46562 Voerde / Germany

Phone: +49 281 - 9 44 00 - 0  
Fax: +49 281 - 9 44 00 - 11  
e-mail: [export@endo-flex.de](mailto:export@endo-flex.de)  
[www.endo-flex.de](http://www.endo-flex.de)



disposable Sclerotherapy Needles **Ø 2.3 mm**,  
for working channels  $\geq \text{Ø } 2.8 \text{ mm}$

needle Ø 0.7 mm



Handle with  
Snap-In-Mechanism



STERILE EO

length	packaging unit	needle Ø 0.7 mm	
		needle length	
		5 mm	
180 cm	5 pieces	NET 25 22-C5	
230 cm	5 pieces	NET 25 22-G5	

disposable Sclerotherapy Needles **Ø 2.3 mm**,  
for working channels  $\geq \text{Ø } 2.8 \text{ mm}$

needle Ø 0.7 mm



needle Ø 0.9 mm



incl. needle safety device ➤



STERILE EO

length	packaging unit	needle Ø 0.7 mm		needle Ø 0.9 mm	
		needle length		needle length	
		4 mm	6 mm	4 mm	6 mm
160 cm	5 pieces		ET 25 22-B6		ET 26 22-B6
180 cm	5 pieces	ET 25 22-C4	ET 25 22-C6		ET 26 22-C6
230 cm	5 pieces	ET 25 22-G4	ET 25 22-G6		ET 26 22-G6
280 cm	5 pieces		ET 25 22-M6		
300 cm	5 pieces		ET 25 22-O6		

disposable Sclerotherapy Needles **Ø 1.8 mm**,  
for working channels  $\geq \text{Ø } 2.0 \text{ mm}$

needle Ø 0.6 mm



incl. needle safety device ➤



STERILE EO

length	packaging unit	needle Ø 0.6 mm	
		needle length	
		4 mm	6 mm
180 cm	5 pieces	ET 25 18-C4	ET 25 18-C6
280 cm	5 pieces		ET 25 18-M6