

TRANSBRONCHIAL ASPIRATION NEEDLES

disposable



 **ENDO-FLEX**[®]

disposable Transbronchial Aspiration Needle \varnothing 1.8 and 2.3 mm

Needle with lateral opening,
 \varnothing 0.7 und 0.9 mm



incl. needle safety device

length	packing unit	needle \varnothing 0.7 mm		needle \varnothing 0.9 mm	
		tube		tube	
		\varnothing 1.8 mm	\varnothing 2.3 mm	\varnothing 1.8 mm	\varnothing 2.3 mm
		needle length 15 mm	needle length 15 mm	needle length 12 mm	needle length 15 mm
120 cm	5 pieces	EB 25 19-A15	EB 25 22-A15	EB 26 18-A12	EB 26 22-A15

Legend:



This instrument is designed and warranted for single use only.



This instrument is sterile packed separately.

ENDO-FLEX GmbH
 Alte Hünxer Str. 115
 46562 Voerde / Germany

Phone: +49 281 - 9 44 00 - 0
 Fax: +49 281 - 9 44 00 - 11
 e-mail: export@endo-flex.de
www.endo-flex.de