


# FIBRIN APPLICATION NEEDLES

re-usable

AUTOCLAVABLE  
134°C  
273°F



 **ENDO-FLEX<sup>®</sup>**

# double lumen probes $\varnothing$ 2.6 mm for Fibrin application, for working channels $\geq \varnothing$ 2.8 mm

double needle  $\varnothing$  0.8 und 0.65 mm,  
needle length: 6 mm



- binding of both components inside of tissue
- simultaneous injection of both components

- colour marked LL connectors:

**RED** = Trombin

**BLUE** = Fibrinogen

length	Order No.
160 cm	T 27 2610-B
180 cm	T 27 2610-C
230 cm	T 27 2610-G

exchangeable needle  $\varnothing$  0.8 mm,  
needle length: 6 mm,  
incl. 1 spare needle



- binding of both components inside of needle tip
- simultaneous injection of both components

- colour marked LL connectors:

**RED** = Trombin

**BLUE** = Fibrinogen

length	Order No.
160 cm	T 27 2620-B
180 cm	T 27 2620-C
230 cm	T 27 2620-G

## Legend:



This instrument is autoclavable by 134°C / 273°F.

This instrument is detachable for cleaning.

ENDO-FLEX GmbH  
Alte Hünxer Str. 115  
46562 Voerde / Germany

Phone: +49 281 - 9 44 00 - 0

Fax: +49 281 - 9 44 00 - 11

e-mail: [export@endo-flex.de](mailto:export@endo-flex.de)

[www.endo-flex.de](http://www.endo-flex.de)