

CYTOLOGY BRUSHES

disposable



ENDO-FLEX[®]

disposable Cytology Brushes \varnothing 1.8 and 2.3 mm

metal ball for closing the tube



length	packing unit	\varnothing 1.8 mm	\varnothing 2.3 mm
		for working channels $\geq \varnothing$ 2.0 mm	for working channels $\geq \varnothing$ 2.8 mm
120 cm	5 pieces	E 42 18-A	
180 cm	5 pieces	E 42 18-C	
230 cm	5 pieces		E 42 22-G

Legend:



This instrument is designed and warranted for single use only.



This instrument is sterile packed separately.

ENDO-FLEX GmbH
Alte Hünxer Str. 115
46562 Voerde / Germany

Phone: +49 281 - 9 44 00 - 0
Fax: +49 281 - 9 44 00 - 11
e-mail: export@endo-flex.de
www.endo-flex.de