

E.R.C.P. CATHETERS

disposable



 **ENDO-FLEX**[®]

E.R.C.P. Catheters \varnothing 1.8 mm for working channels $\geq \varnothing$ 2.0 mm, accepts guide wire .035 inch

radiopaque tip:



type	standard tip	tapered tip	hour glass metal tip
accepts guide wire .035 inch	E 43 18-GW	E 44 18-GW	E 48 18-GW

recommended equipment

Appropriate guide wires can be found in chapter 06.



Legend:



This instrument is designed and warranted for single use only.



This instrument is sterile packed separately.

ENDO-FLEX GmbH
Alte Hünxer Str. 115
46562 Voerde / Germany

Phone: +49 281 - 9 44 00 - 0
Fax: +49 281 - 9 44 00 - 11
e-mail: export@endo-flex.de
www.endo-flex.de