


E.R.C.P. Catheters **Ø 1.8 mm** for working channels $\geq \text{Ø } 2.0 \text{ mm}$,
not to be used with guide wire

needle tip 2 mm




type	3 black markers	coloured marker	connector
with mandrin	49 18	49 181	

E.R.C.P. Catheters **Ø 1.6 mm** for working channels $\geq \text{Ø } 2.0 \text{ mm}$,
 accepts guide wire .025 inch

tapered tip



execution	3 black markers	coloured markers	connector
with mandrin	44 16	44 161	

recommended equipment

Appropriate guide wires can be found in chapter 06.



© ENDO-FLEX - Technical specifications are subject to change

Legend:



This instrument is autoclavable by 134°C / 273°F.

ENDO-FLEX GmbH
 Alte Hünxer Str. 115
 46562 Voerde / Germany
 Phone: +49 281 - 9 44 00 - 0
 Fax: +49 281 - 9 44 00 - 11
 e-mail: export@endo-flex.de
 www.endo-flex.de


E.R.C.P. CATHETERS

re-usable




E.R.C.P. Catheters Ø 1.8 mm,
accepts guide wire .035 inch

standard tip



type	3 black markers	coloured markers	connector
with mandrin for working channels ≥ Ø 2.0 mm	43 18	43 181	
with mandrin / proximal Ø 2.45 mm for working channels ≥ Ø 2.8 mm	43 22-18FD	43 221-18FD	


tapered tip



type	3 black markers	coloured markers	connector
with mandrin for working channels ≥ Ø 2.0 mm	44 18	44 181	
with mandrin / proximal Ø 2.45 mm for working channels ≥ Ø 2.8 mm	44 22-18FD	44 221-18FD	


E.R.C.P. Catheters Ø 1.8 mm,
accepts guide wire .035 inch

conical
metal tip



type	3 black markers	coloured markers	connector
with mandrin for working channels ≥ Ø 2.0 mm	45 18	45 181	
with mandrin / proximal Ø 2.45 mm for working channels ≥ Ø 2.8 mm	45 22-18FD	45 221-18FD	


tapered
metal tip



type	3 black markers	coloured markers	connector
with mandrin for working channels ≥ Ø 2.0 mm	46 18	46 181	
with mandrin / proximal Ø 2.45 mm for working channels ≥ Ø 2.8 mm	46 22-18FD	46 221-18FD	


E.R.C.P. Catheters Ø 1.8 mm,
accepts guide wire .035 inch

round
metal tip



type	3 black markers	coloured markers	connector
with mandrin for working channels ≥ Ø 2.0 mm	47 18	47 181	
with mandrin / proximal Ø 2.45 mm for working channels ≥ Ø 2.8 mm	47 22-18FD	47 221-18FD	

hour glass
metal tip



type	3 black markers	coloured markers	connector
with mandrin for working channels ≥ Ø 2.0 mm	48 18	48 181	
with mandrin / proximal Ø 2.45 mm for working channels ≥ Ø 2.8 mm	48 22-18FD	48 221-18FD	