

COAGULATION PROBES

re-usable

AUTOCLAVABLE
134°C
273°F



 **ENDO-FLEX**[®]

Coagulation Probe Ø 1.8 and 2.3 mm

Tip with 1 opening,
handle with flushing port



length	Ø 1.8 mm for working channels ≥ Ø 2.0 mm	Ø 2.3 mm for working channels ≥ Ø 2.8 mm
120 cm	901 80 18-A	901 80 22-A
160 cm	901 80 18-B	901 80 22-B
180 cm	901 80 18-C	901 80 22-C
230 cm	901 80 18-G	901 80 22-G

recommended equipment

active cord:



HF generator	Order No.
ERBE E- and T-Series	64 0300
ERBE ICC	64 0500
OLYMPUS HF 120/300	64 0500



Legend:



This instrument is autoclavable by 134°C / 273°F.



This Instrument possesses a flushing port.



This instrument is designed to be used with HF current.

ENDO-FLEX GmbH
 Alte Hünxer Str. 115
 46562 Voerde / Germany
 Phone: +49 281 - 9 44 00 - 0
 Fax: +49 281 - 9 44 00 - 11
 e-mail: export@endo-flex.de
www.endo-flex.de