

disposable Needle-Knife Sphincterotomes, fixed handle



pointed type,
stop at 3 mm



accepts guide wire
.035 inch
pointed type,
stop at 6 mm

single lumen

Ø 1.8 mm
for working channels ≥ Ø 2.0 mm



OE 110 18N3

double lumen

Ø 2.2 mm distal Ø 1.8 mm
for working channels ≥ Ø 2.8 mm



OE 110 18N-DL



blunt type,
stop at 3 mm

single lumen

Ø 1.8 mm
for working channels ≥ Ø 2.0 mm
OE 110 18SP3

recommended equipment


active cords:



HF generator	Order No.
ERBE E- and T-Series	64 0300
ERBE ICC	64 0500
OLYMPUS HF 120/300	64 0500



Appropriate guide wires can be found in chapter 06.



© ENDO-FLEX - Technical specifications are subject to change

Legend:

- This instrument is designed and warranted for single use only.
- This instrument is designed to be used with HF current.
- This instrument is sterile packed separately.
- This instrument is autoclavable by 134°C / 273°F.

ENDO-FLEX GmbH
Alte Hünxer Str. 115
46562 Voerde / Germany
Phone: +49 281 - 9 44 00 - 0
Fax: +49 281 - 9 44 00 - 11
e-mail: export@endo-flex.de
www.endo-flex.de

SPHINCTEROTOMES
disposable




disposable **non-guide-wire** Sphincterotomes **Ø 1.8 mm**, fixed handle,
for working channels $\geq \text{Ø } 2.0 \text{ mm}$



tapered tip

single lumen

cutting length	20 mm	30 mm
	OE 105 1821	OE 105 1831



precut type

single lumen



precut type
with tapered tip

single lumen

cutting length	20 mm	30 mm	20 mm	30 mm
	OE 103 1821	OE 103 1831	OE 105 1821PC	OE 105 1831PC


disposable **guide wire** Sphincterotomes, **Ø 2.2 mm** distal **Ø 1.8 mm**,
fixed handle, for working channels $\geq 2.8 \text{ mm } \text{Ø}$



standard tip

double lumen

accepts guide wire .035 inch	double lumen			double lumen		
cutting length	20 mm	25 mm	30 mm	20 mm	25 mm	30 mm
	OE 100 2220DL	OE 100 2225DL	OE 100 2230DL	OE 102 2220DL	OE 102 2225DL	OE 102 2230DL



precut type

double lumen



tapered tip

double lumen

accepts guide wire .035 inch	double lumen			double lumen		
cutting length	20 mm	25 mm	30 mm	20 mm	25 mm	30 mm
	OE 104 2220DL	OE 104 2225DL	OE 104 2230DL	OE 104 2220DLC	OE 104 2225DLC	OE 104 2230DLC




„SmartCut“
partially insulated,
tapered tip


double lumen

recommended equipment

Appropriate guide wires can be found in chapter 06.



disposable **guide wire** Sphincterotomes, **Ø 2.2 mm** distal **Ø 1.8 mm**,
fixed handle, for working channels $\geq 2.8 \text{ mm } \text{Ø}$



precut type

accepts guide wire
.035 inch

triple lumen

cutting length	20 mm	25 mm	30 mm
	OE 102 2220TRL	OE 102 2225TRL	OE 102 2230TRL



tapered tip

triple lumen

accepts guide wire .035 inch	triple lumen			triple lumen		
cutting length	20 mm	25 mm	30 mm	20 mm	25 mm	30 mm
	OE 104 2220TRL	OE 104 2225TRL	OE 104 2230TRL	OE 104 2220TRLC	OE 104 2225TRLC	OE 104 2230TRLC



„SmartCut“
partially insulated,
tapered tip

triple lumen

recommended equipment

Appropriate guide wires can be found in chapter 06.

