

# BILIARY DILATION CATHETERS

re-usable


AUTOCLAVABLE  
134°C  
273°F



 **ENDO-FLEX**<sup>®</sup>

# Biliary Dilation Catheters

- catheter length 180 cm
- radiopaque tube
- distal radiopaque markers



catheter tip	Fr. 4.0	Fr. 5.0	Fr. 5.0	Fr. 5.0	Fr. 5.0	Fr. 7.0	Fr. 7.0
catheter Ø	Fr. 6.0	Fr. 7.0	Fr. 8.5	Fr. 9.0	Fr. 10.0	Fr. 10.0	Fr. 12.0
accepts guide wire	.035 inch	.035 inch	.035 inch	.035 inch	.035 inch	.035 inch	.035 inch
Order No.	200 00064	200 00075	200 00085	200 00095	200 00105	200 00107	200 00127
colour mark	violet	yellow	orange	white	red	red	blue

## recommended equipment

Appropriate guide wires can be found in chapter 06.



### Legend:



This Instrument possesses a flushing port.



This instrument is autoclavable by 134°C / 273°F.

ENDO-FLEX GmbH  
Alte Hünxer Str. 115  
46562 Voerde / Germany

Phone: +49 281 - 9 44 00 - 0  
Fax: +49 281 - 9 44 00 - 11  
e-mail: [export@endo-flex.de](mailto:export@endo-flex.de)  
[www.endo-flex.de](http://www.endo-flex.de)