






# P T C D



© ENDO-FLEX - Technical specifications are subject to change

## Legend:

-  This instrument is designed and warranted for single use only.
-  This instrument is sterile packed separately.
-  This instrument is autoclavable by 134°C / 273°F.
-  This Instrument possesses a flushing port.

ENDO-FLEX GmbH  
Alte Hünxer Str. 115  
46562 Voerde / Germany  
Phone: +49 281 - 9 44 00 - 0  
Fax: +49 281 - 9 44 00 - 11  
e-mail: [export@endo-flex.de](mailto:export@endo-flex.de)  
[www.endo-flex.de](http://www.endo-flex.de)

 **ENDO-FLEX**®

## disposable PTC-D Set Fr. 7.0 / consisting of 4 parts

PTC needle



ring-drainage-catheter



guide wire



PTC exchange guide wire



Order No.	name
280 170	PTCD-Set Fr. 7.0 / consisting of 4 parts

single components:

280 171	PTC needle 19 G incl. catheter Fr. 5.0 / length 30 cm
280 172	teflon coated stainless steel guide wire .035 inch / length 100 cm
280 173	PTC exchanging guide wire .035 inch / 7.5 cm flexible part / length 82.5 cm
280 174	ring-drainage-catheter Fr. 7.0 / 6 sideports / length 50 cm



STERILE EO

## disposable PTC-D Set Fr. 8.4 / consisting of 6 parts


PTC needle



dilator



guide wire



ring-drainage-catheter



PTC exchange guide wire



drainage bag connector



Order No.	name
280 180	PTCD-Set Fr. 8.4 / consisting of 6 parts

single components:


280 182	PTC needle 19 G incl. catheter Fr. 5.0 / length 25 cm
280 183	teflon coated guide wire .035 inch / 3 mm J-curved / length 80 cm
280 184	PTC exchanging guide wire .035 inch / 7.5 cm flexible part / length 82.5 cm
280 185	dilator Fr. 8.0 / length 20 cm
280 186	ring-drainage-catheter Fr. 8.4 / 32 sideports / length 50 cm
280 187	drainage bag connector / funnel-male / length 30 cm



STERILE EO

## re-usable PTCD Dilation Catheters

- catheter length 40 cm
- radiopaque tube
- distal radiopaque markers



catheter tip	Fr. 4.0	Fr. 5.0	Fr. 5.0	Fr. 5.0	Fr. 5.0	Fr. 7.0	Fr. 7.0
catheter Ø	Fr. 6.0	Fr. 7.0	Fr. 8.5	Fr. 9.0	Fr. 10.0	Fr. 10.0	Fr. 12.0
accepts guide wire	.035 inch	.035 inch	.035 inch	.035 inch	.035 inch	.035 inch	.035 inch
Order No.	200 00064-40	200 00075-40	200 00085-40	200 00095-40	200 00105-40	200 00107-40	200 00127-40
colour mark	violet	yellow	orange	white	red	red	blue

## recommended equipment

Appropriate guide wires can be found in chapter 06.

