

## kink resistant Memory Basket $\varnothing$ 2.3 mm & $\varnothing$ 1.8 mm

basket made of gold plated Nitinol,  
radiopaque tube,  
incl. handle:



4 wires



	dormia-type			basket type		
	basket size (in stretched condition)			basket size (in stretched condition)		
$\varnothing$ 2.3 mm for working channels $\geq \varnothing$ 2.8 mm	40 mm	50 mm	60 mm	40 mm	50 mm	60 mm
Nitinol wire	146 224	146 225	146 226	141 224	141 225	141 226

	dormia-type			basket type		
	basket size (in stretched condition)			basket size (in stretched condition)		
$\varnothing$ 1.8 mm for working channels $\geq \varnothing$ 2.0 mm	40 mm	50 mm	60 mm	40 mm	50 mm	60 mm
Nitinol wire	146 184	146 185	146 186	141 184	141 185	141 186

6 wires



	dormia-type			dormia-type		
	basket size (in stretched condition)			basket size (in stretched condition)		
$\varnothing$ 2.3 mm for working channels $\geq \varnothing$ 2.8 mm	40 mm	50 mm	60 mm	40 mm	50 mm	60 mm
Nitinol wire	147 224	147 225	147 226	142 224	142 225	142 226

	dormia-type			dormia-type		
	basket size (in stretched condition)			basket size (in stretched condition)		
$\varnothing$ 1.8 mm for working channels $\geq \varnothing$ 2.8 mm	40 mm	50 mm	60 mm	40 mm	50 mm	60 mm
Nitinol wire	147 184	147 185	147 186	142 184	142 185	142 186

NOTE: Kink resistant Memory Baskets are not suitable for lithotripsy!

## kink resistant Memory Basket $\varnothing$ 2.6 mm, accepts guide wire .035 inch, for working channels $\geq \varnothing$ 2.8 mm

basket made of gold plated Nitinol,  
radiopaque tube,  
incl. handle:



4 wires



	dormia-type			basket type		
	basket size (in stretched condition)			basket size (in stretched condition)		
	40 mm	50 mm	60 mm	40 mm	50 mm	60 mm
Nitinol wire	146 264GW	146 265GW	146 266GW	141 264GW	141 265GW	141 266GW

NOTE: Kink resistant Memory Baskets are not suitable for lithotripsy!

### recommended equipment

Appropriate guide wires can be found in chapter 06.



#### Legend:

- This instrument is autoclavable by 134°C / 273°F.
- This instrument is rotatable.
- This instrument is detachable for cleaning.

ENDO-FLEX GmbH  
Alte Hünxer Str. 115  
46562 Voerde / Germany  
Phone: +49 281 - 9 44 00 - 0  
Fax: +49 281 - 9 44 00 - 11  
e-mail: export@endo-flex.de  
www.endo-flex.de

# EXTRACTION BASKETS

re-usable




Extraction Baskets  $\varnothing$  2.3 mm for working channels  $\geq \varnothing$  2.8 mm

basket made of stainless steel, incl. handle:




4 wires



	dormia-type				basket-type			
	basket size (in stretched condition)				basket size (in stretched condition)			
	40 mm	50 mm	60 mm	70 mm	40 mm	50 mm	60 mm	70 mm
monofilament wire	150 224	150 225	150 226	150 227	151 224	151 225	151 226	151 227
braided wire	160 224	160 225	160 226	160 227	161 224	161 225	161 226	161 227

6 wires



	dormia-type				basket-type			
	basket size (in stretched condition)				basket size (in stretched condition)			
	40 mm	50 mm	60 mm	70 mm	40 mm	50 mm	60 mm	70 mm
monofilament wire	152 224	152 225	152 226	152 227	153 224	153 225	153 226	153 227
braided wire	162 224	162 225	162 226	162 227	163 224	163 225	163 226	163 227

Extraction Baskets  $\varnothing$  2.3 mm, accepts guide wire .035 inch, for working channels  $\geq \varnothing$  2.8 mm

basket made of stainless steel, incl. handle:



4 wires



	dormia-type				basket-type			
	basket size (in stretched condition)				basket size (in stretched condition)			
	40 mm	50 mm	60 mm	70 mm	40 mm	50 mm	60 mm	70 mm
monofilament wire	150 224GW	150 225GW	150 226GW	150 227GW	151 224GW	151 225GW	151 226GW	151 227GW
braided wire	160 224GW	160 225GW	160 226GW	160 227GW	161 224GW	161 225GW	161 226GW	161 227GW

6 wires




	dormia-type				basket-type			
	basket size (in stretched condition)				basket size (in stretched condition)			
	40 mm	50 mm	60 mm	70 mm	40 mm	50 mm	60 mm	70 mm
monofilament wire	152 224GW	152 225GW	152 226GW	152 227GW	153 224GW	153 225GW	153 226GW	153 227GW

Extraction Baskets  $\varnothing$  1.8 mm for working channels  $\geq \varnothing$  2.0 mm

basket made of stainless steel, incl. handle:




4 wires



	dormia-type			basket-type		
	basket size (in stretched condition)			basket size (in stretched condition)		
	30 mm	40 mm	50 mm	30 mm	40 mm	50 mm
monofilament wire	150 183	150 184	150 185	151 183	151 184	151 185
braided wire	160 183	160 184	160 185	161 183	161 184	161 185

6 wires



	dormia-type			basket-type		
	basket size (in stretched condition)			basket size (in stretched condition)		
	30 mm	40 mm	50 mm	30 mm	40 mm	50 mm
monofilament wire	152 183	152 184	152 185	153 183	153 184	153 185

recommended equipment

Appropriate guide wires can be found in chapter 06.

