

# LITHOTRIPTOR BASKETS

disposable



 **ENDO-FLEX**<sup>®</sup>

# Lithotripter Baskets $\varnothing$ 2.6 mm for working channel $\geq \varnothing$ 2.8 mm

- tube with lateral contrast injection port
- braided traction cable
- total length 4 m



STERILE EO

braided wire	4 wires - basket type		6 wires - dormia type	
	basket size (in stretched condition)		basket size (in stretched condition)	
	50 mm	60 mm	50 mm	60 mm
	E 122 265N	E 122 266N	E 124 265N	E 124 266N

to be used with:

- mechanical Lithotripter Set Type I  
(see following pages)



- mechanical Lithotripter Set Type III  
(see following pages)



## Legend:



This instrument is designed and warranted for single use only.



This instrument is sterile packed separately.



This instrument is autoclavable by 134°C / 273°F.

ENDO-FLEX GmbH  
 Alte Hünxer Str. 115  
 46562 Voerde / Germany  
 Phone: +49 281 - 9 44 00 - 0  
 Fax: +49 281 - 9 44 00 - 11  
 e-mail: export@endo-flex.de  
 www.endo-flex.de