

Lithotripter Baskets \varnothing 2.6 mm for working channel $\geq \varnothing$ 2.8 mm
accepts guide wire .035 inch

- tube with lateral contrast injection port
- braided traction cable
- total length 4 m



4 wires
basket type



6 wires
dormia type



4 wires	basket type			
	basket size (in stretched condition)			
	50 mm	60 mm	70 mm	80 mm
braided wire	122 265N-GW	122 266N-GW	122 267N-GW	122 268N-GW

6 wires	dormia type			
	basket size (in stretched condition)			
	50 mm	60 mm	70 mm	80 mm
braided wire	124 265N-GW	124 266N-GW	124 267N-GW	124 268N-GW

recommended equipment

	handle	stainless steel spiral
Order No.	120 000N	for working channels $\geq \varnothing$ 2.8 mm: 120 261 (\varnothing 2.6 mm / length 195 cm)
Order No.	120 000N	for working channels $\geq \varnothing$ 3.7 mm: 120 331 (\varnothing 3.4 mm / length 195 cm)

Appropriate guide wires can be found in chapter 06.



Legend:



This instrument is autoclavable by 134°C / 273°F.



This instrument is detachable for cleaning.

ENDO-FLEX GmbH
Alte Hünxer Str. 115
46562 Voerde / Germany

Phone: +49 281 - 9 44 00 - 0
Fax: +49 281 - 9 44 00 - 11
e-mail: export@endo-flex.de
www.endo-flex.de

MECHANICAL LITHOTRIPTORS
re-usable



© ENDO-FLEX - Technical specifications are subject to change

Mechanical Lithotripter Set Type I \varnothing 2.6 mm

handle with coaxial fixation, stainless steel spiral, 2 baskets

For methodical lithotripsy „Through-The-Scope“.



baskets made of braided wire, with 4 m traction cable

tube with lateral contrast injection port

stainless steel spiral \varnothing 2.6 mm

handle with coaxial fixation


Set Type I	basket size (in stretched condition)		
	50 mm	60 mm	70 mm
Order No.	120 265N	120 266N	120 267N

Replacement

	handle	stainless steel spiral
Order No.	120 000N	for working channels $\geq \varnothing$ 2.8 mm: 120 261 (\varnothing 2.6 mm / length 195 cm)

Lithotripter Baskets \varnothing 2.6 mm for working channel $\geq \varnothing$ 2.8 mm

- tube with lateral contrast injection port and radiopaque tip
- braided traction cable
- total length 4 m



4 wires basket type

6 wires basket type

6 wires dormia type

4 wires	basket type		
	basket size (in stretched condition)		
	50 mm	60 mm	70 mm
braided wire	122 265N	122 266N	122 267N

6 wires	dormia type			basket type		
	basket size (in stretched condition)			basket size (in stretched condition)		
	50 mm	60 mm	70 mm	50 mm	60 mm	70 mm
braided wire	124 265N	124 266N	124 267N	128 265N	128 266N	128 267N

Lithotripter Baskets \varnothing 3.4 mm for working channel $\geq \varnothing$ 3.7 mm

4 wires	basket type		
	basket size (in stretched condition)		
	50 mm	60 mm	70 mm
monofilament wire	121 335N	121 336N	121 337N
braided wire	122 335N	122 336N	122 337N

6 wires	dormia type			basket type		
	basket size (in stretched condition)			basket size (in stretched condition)		
	50 mm	60 mm	70 mm	50 mm	60 mm	70 mm
braided wire	124 335N	124 336N	124 337N	128 335N	128 336N	128 337N

recommended equipment

	handle	stainless steel spiral
Order No.	120 000N	for working channels $\geq \varnothing$ 2.8 mm: 120 261 (\varnothing 2.6 mm / length 195 cm)
Order No.	120 000N	for working channels $\geq \varnothing$ 3.7 mm: 120 331 (\varnothing 3.4 mm / length 195 cm)