

MECHANICAL LITHOTRIPTOR SET TYPE III re-usable

AUTOCLAVABLE
134°C
273°F



 **ENDO-FLEX**®

Mechanical Lithotripter Set Type III in storage case

Combination of methodical and emergency lithotripter in one device.

Easy handling by unique spiral and traction cable fixation.

Compatible to all common ENDO-FLEX baskets.



Order No.	components
120 610	1 x storage case 1 x handle 2 x stainless steel spiral Ø 2.6 mm / length 195 cm 1 x stainless steel spiral Ø 4.0 mm / length 100 cm 1 x basket, Ø 2.6 mm, 4 braided wires, basket type, basket size 50 mm 1 x basket, Ø 2.6 mm, 4 braided wires, basket type, basket size 60 mm 1 x basket, Ø 2.6 mm, 4 braided wires, basket type, basket size 70 mm


Single Components



handle



stainless steel spiral



basket



Order No.	name
120 602	stainless steel spiral Ø 2.6 mm / length 195 cm
120 603	stainless steel spiral Ø 3.4 mm / length 195 cm
120 604	stainless steel spiral Ø 4.0 mm / length 100 cm
120 611	handle incl. 2 x stainless steel spirals (1 x Order No. 120 602 and 1 x Order No. 120 604)
122 265N	basket, Ø 2.6 mm, 4 braided wires, basket type, basket size 50 mm
122 266N	basket, Ø 2.6 mm, 4 braided wires, basket type, basket size 60 mm
122 267N	basket, Ø 2.6 mm, 4 braided wires, basket type, basket size 70 mm

Legend:



This instrument is autoclavable by 134°C / 273°F.



This instrument is detachable for cleaning.

ENDO-FLEX GmbH
 Alte Hünxer Str. 115
 46562 Voerde / Germany
 Phone: +49 281 - 9 44 00 - 0
 Fax: +49 281 - 9 44 00 - 11
 e-mail: export@endo-flex.de
 www.endo-flex.de