

PROGRESSIVE BALLOON DILATOR


disposable



 **ENDO-FLEX[®]**

Progressive Balloon Dilator

- catheter length 230 cm
- radiopaque markers
- 3 progressive stages



Two types of application combined in one balloon:

1. TTS application via stylet.
2. OTW application over a guide wire .035 inch after removing the stylet.


Order No.	Ø catheter	balloon length	fitting for working channel	Ø balloon		
				3 atm	4.5 atm	6 atm
341 15PRO	7 Fr.	55 mm	2.8 mm	15 mm	16.5 mm	18 mm
341 18PRO	7 Fr.	55 mm	2.8 mm	18 mm	19 mm	20 mm

recommended equipment

Inflation / Deflation Device incl. PSI-Manometer 0-440 PSI



Order No.	340 800 INFL
-----------	--------------



Appropriate guide wires can be found in chapter 06.

Legend:



This instrument is designed and warranted for single use only.



This instrument is sterile packed separately.

ENDO-FLEX GmbH
 Alte Hünxer Str. 115
 46562 Voerde / Germany
 Phone: +49 281 - 9 44 00 - 0
 Fax: +49 281 - 9 44 00 - 11
 e-mail: export@endo-flex.de
 www.endo-flex.de