

ESOPHAGEAL DILATOR SETS

re-usable



 **ENDO-FLEX[®]**

Esophageal Dilator Sets

- flexible, radiopaque plastic
- accepts guide wire
- coloured marks for exact measurement


Set 1:

Order No.	components
310 070-1	7 dilators: Ø 5-7-9-11-13-14-15 mm / length 70 cm 1 guide wire 1 cleaning brush 25 disposable flushing adapters 1 storage case



Set 2:

Order No.	components
310 070-2	16 dilators: Ø 5 - 20 mm / length 70 cm 1 guide wire 1 cleaning brush 25 disposable flushing adapters 1 storage case



Single Components

Dilators

Order No.	parameters
310 070-05	dilator Fr. 15 / Ø 5 mm / length 70 cm
310 070-06	dilator Fr. 18 / Ø 6 mm / length 70 cm
310 070-07	dilator Fr. 21 / Ø 7 mm / length 70 cm
310 070-08	dilator Fr. 24 / Ø 8 mm / length 70 cm
310 070-09	dilator Fr. 27 / Ø 9 mm / length 70 cm
310 070-10	dilator Fr. 30 / Ø 10 mm / length 70 cm
310 070-11	dilator Fr. 33 / Ø 11 mm / length 70 cm
310 070-12	dilator Fr. 36 / Ø 12 mm / length 70 cm

Order No.	parameters
310 070-13	dilator Fr. 39 / Ø 13 mm / length 70 cm
310 070-14	dilator Fr. 42 / Ø 14 mm / length 70 cm
310 070-15	dilator Fr. 45 / Ø 15 mm / length 70 cm
310 070-16	dilator Fr. 48 / Ø 16 mm / length 70 cm
310 070-17	dilator Fr. 51 / Ø 17 mm / length 70 cm
310 070-18	dilator Fr. 54 / Ø 18 mm / length 70 cm
310 070-19	dilator Fr. 57 / Ø 19 mm / length 70 cm
310 070-20	dilator Fr. 60 / Ø 20 mm / length 70 cm

Guide Wire Ø 1.8 mm with 5 cm flexible tip

Order No.	characteristic
311 100	non marked
311 100M	marked



Cleaning Brush

Order No.	packing unit
310 070-RB	3 pieces



Flushing Adapter

Order No.	packing unit
300 070-RA	25 pieces



Legend:



This instrument is designed and warranted for single use only.



This instrument is autoclavable by 134°C / 273°F.

ENDO-FLEX GmbH
 Alte Hünxer Str. 115
 46562 Voerde / Germany
 Phone: +49 281 - 9 44 00 - 0
 Fax: +49 281 - 9 44 00 - 11
 e-mail: export@endo-flex.de
 www.endo-flex.de