

PNEUMATIC DILATORS

re-usable




 **ENDO-FLEX**[®]

Pneumatic Dilator by Witzel for treatment of stenosis in esophagus, stomach and duodenum

- max. endoscope Ø 11.2 mm
- dilator up to max. Ø 20 mm
- incl. inflation device up to 300 mm Hg



dilator length	Order No.
20 cm	340 650




Pneumatic Achalasia Dilator by Witzel for treatment of achalasia

- max. endoscope Ø 11.2 mm
- dilator up to max. Ø 40 mm
- incl. inflation device up to 300 mm Hg



dilator length	Order No.
16 cm	340 660



Legend:



This instrument is sterilizable by gamma-ray.

ENDO-FLEX GmbH
Alte Hünxer Str. 115
46562 Voerde / Germany
Phone: +49 281 - 9 44 00 - 0
Fax: +49 281 - 9 44 00 - 11
e-mail: export@endo-flex.de
www.endo-flex.de