


self-expanding Colonic Stents

- incl. pre-loaded application set
- distal placement
- radiopaque markers




uncoated



length	stent diameter	tulip diameter	type	Order No.
100 mm	20 mm	24 mm	uncoated	COL-0-20-100
120 mm	20 mm	24 mm	uncoated	COL-0-20-120
80 mm	22 mm	26 mm	uncoated	COL-0-22- 80
100 mm	22 mm	26 mm	uncoated	COL-0-22-100

Application Set (not available separately)




for type	diameter	length	accepts guide wire
uncoated	Fr. 11.0	180 cm	.035 inch

recommended equipment

Appropriate guide wires can be found in chapter 06.



NITINOL-STENTS disposable




© ENDO-FLEX - Technical specifications are subject to change

Legend:



This instrument is designed and warranted for single use only.




This instrument is sterile packed separately.

ENDO-FLEX GmbH
 Alte Hünxer Str. 115
 46562 Voerde / Germany
 Phone: +49 281 - 9 44 00 - 0
 Fax: +49 281 - 9 44 00 - 11
 e-mail: export@endo-flex.de
 www.endo-flex.de



self-expanding Biliary Stents


- incl. pre-loaded application set
- trans-papillary (proximal) placement
- radiopaque markers



uncoated



length	stent diameter	type	Order No.
60 mm	10 mm	uncoated	BIL-0-10- 60-RP
80 mm	10 mm	uncoated	BIL-0-10- 80-RP
100 mm	10 mm	uncoated	BIL-0-10-100-RP




silicone coated



length	stent diameter	type	Order No.
60 mm	10 mm	coated	BIL-1-10-60-RP
80 mm	10 mm	coated	BIL-1-10-80-RP


Application Set (not available separately)



for type	diameter	length	accepts guide wire
uncoated	Fr. 9.0	180 cm	.035 inch
coated	Fr. 10.0	180 cm	.035 inch

self-expanding Esophageal Stents

- incl. pre-loaded application set
- radiopaque markers




silicone coated



	length	stent diameter	tulip diameter	type	Order No.
distal placement	80 mm	20 mm	24 mm	partially coated	ESO-1-20- 80
	100 mm	20 mm	24 mm	partially coated	ESO-1-20-100
	120 mm	20 mm	24 mm	partially coated	ESO-1-20-120
	140 mm	20 mm	24 mm	partially coated	ESO-1-20-140
distal placement	80 mm	24 mm	28 mm	partially coated	ESO-1-24- 80
	100 mm	24 mm	28 mm	partially coated	ESO-1-24-100
	120 mm	24 mm	28 mm	partially coated	ESO-1-24-120
	140 mm	24 mm	28 mm	partially coated	ESO-1-24-140
proximal placement	80 mm	20 mm	24 mm	partially coated	ESO-1-20- 80-RP
	100 mm	20 mm	24 mm	partially coated	ESO-1-20-100-RP
	120 mm	20 mm	24 mm	partially coated	ESO-1-20-120-RP
	140 mm	20 mm	24 mm	partially coated	ESO-1-20-140-RP
proximal placement	80 mm	24 mm	28 mm	partially coated	ESO-1-24- 80-RP
	100 mm	24 mm	28 mm	partially coated	ESO-1-24-100-RP
	120 mm	24 mm	28 mm	partially coated	ESO-1-24-120-RP
	140 mm	24 mm	28 mm	partially coated	ESO-1-24-140-RP


Application Set (not available separately)



for type	diameter	length	accepts guide wire
coated	Fr. 24.0	70 cm	.035 inch

self-expanding Duodenal Stents

- incl. pre-loaded application set
- distal placement
- radiopaque markers



uncoated



length	stent diameter	tulip diameter	type	Order No.
80 mm	20 mm	24 mm	uncoated	DUO-0-20- 80
100 mm	20 mm	24 mm	uncoated	DUO-0-20-100

Application Set (not available separately)



for type	diameter	length	accepts guide wire
uncoated	Fr. 10.0	180 cm	.035 inch