

COLON DECOMPRESSION PROBES

disposable



 **ENDO-FLEX[®]**

disposable Colon Decompression Probes

- accepts guide wire \leq .052 inch
- 6 distal sideports
- universal connector
- radiopaque tube



tube	length	Order No.
Fr. 14.0	175 cm	291 140
Fr. 16.0	175 cm	291 160
Fr. 18.0	175 cm	291 180

recommended equipment

teflon coated stainless steel Guide Wire



Order No.	diameter	length	flexible part
210 52400	.052 inch	400 cm	3 cm

Legend:



This instrument is designed and warranted for single use only.



This instrument is sterile packed separately.

ENDO-FLEX GmbH
Alte Hünxer Str. 115
46562 Voerde / Germany

Phone: +49 281 - 9 44 00 - 0
Fax: +49 281 - 9 44 00 - 11
e-mail: export@endo-flex.de
www.endo-flex.de