


ENDO-NEWS

Progressive Over-The-Wire Balloon Dilator - 3 progressive stages -

- applicable for esophageal or colonic treatments
- 7 Fr. sheath diameter only
- fitting to working channel 2.8 mm
- for guide wire .035 inch
- balloon length 55 mm
- catheter length 230 cm

NEW! Two types of application combined in one balloon:

1. TTS application via the inserted stylet.
2. OTW application over a guide wire .035" after removing the stylet.



Order No.	catheter-Ø	catheter length	balloon-Ø (under pressure of)			balloon length
34106PRO	7 Fr.	230 cm	6 mm (3.0 atm)	7 mm (6.0 atm)	8 mm (10.0 atm)	55 mm
34108PRO	7 Fr.	230 cm	8 mm (3.0 atm)	9 mm (5.5 atm)	10 mm (8.0 atm)	55 mm
34110PRO	7 Fr.	230 cm	10 mm (3.0 atm)	11 mm (5.0 atm)	12 mm (7.0 atm)	55 mm
34112PRO	7 Fr.	230 cm	12 mm (3.0 atm)	13.5 mm (4.5 atm)	15 mm (7.5 atm)	55 mm
34115PRO	7 Fr.	230 cm	15 mm (3.0 atm)	16.5 mm (4.5 atm)	18 mm (6.5 atm)	55 mm
34118PRO	7 Fr.	230 cm	18 mm (3.0 atm)	19 mm (4.5 atm)	20 mm (5.5 atm)	55 mm

For more information about the products please contact your distributor.

04 / 2012



ENDO-FLEX®

ENDO-FLEX GmbH • Alte Hünxer Straße 115 • 46562 Voerde • Germany

☎ +49 (0)281 / 9 44 00 - 0
Fax. +49 (0)281 / 9 44 00 - 11
e-mail: export@endo-flex.de
<http://www.endo-flex.de>