


ENDO-NEWS

Disposable Biliary Dilation Catheter




- radiopaque teflon tube
- accepts guide wire .035"
- radiopaque markers
- available in two lengths
- individually packed
- packing unit : 5 pieces

CE0483



STERILE EO

Order No.	Working Length	Usage	Catheter Diam.	accepts guide wire
E20000064	180 cm	endoscopic	Fr. 6 - Fr. 4	.035 "
E20000064-40	40 cm	PTCD		
E20000075	180 cm	endoscopic	Fr. 7 - Fr. 5	.035 "
E20000075-40	40 cm	PTCD		
E20000085	180 cm	endoscopic	Fr. 8 - Fr. 5	.035 "
E20000085-40	40 cm	PTCD		
E20000095	180 cm	endoscopic	Fr. 9 - Fr. 5	.035 "
E20000095-40	40 cm	PTCD		
E20000105	180 cm	endoscopic	Fr. 10 - Fr. 5	.035 "
E20000105-40	40 cm	PTCD		
E20000107	180 cm	endoscopic	Fr. 10 - Fr. 7	.035 "
E20000107-40	40 cm	PTCD		
E20000127	180 cm	endoscopic	Fr. 12 - Fr. 7	.035 "
E20000127-40	40 cm	PTCD		



For more information about the products please contact your distributor.

02 / 2011



ENDO-FLEX®

ENDO-FLEX GmbH • Alte Hünxer Straße 115 • 46562 Voerde • Germany

☎ +49 (0)281 / 9 44 00 - 0
Fax. +49 (0)281 / 9 44 00 - 11
e-mail: export@endo-flex.de
<http://www.endo-flex.de>