



ENDO-NEWS

Disposable Transbronchial Aspiration Needle



- outer metal sheath for enhanced punction force
- exposed needle length 12 mm
- needle diameter 0.8 mm / 21 G
- lateral opening
- Click-Mechanism for Single-Hand-Control
- single packed sterile
- for single-use only



STERILE EO



CE0483

Order No.	length	packing unit
MEB2618-A12	120 cm	5 pieces

For more information about the products please contact your distributor.

09 / 2009



ENDO-FLEX®

ENDO-FLEX GmbH • Alte Hünxer Straße 115 • 46562 Voerde • Germany

☎ +49 (0)281 / 9 44 00 - 0
Fax. +49 (0)281 / 9 44 00 - 11
e-mail: export@endo-flex.de
<http://www.endo-flex.de>