


ENDO-NEWS

LATEX-FREE Biliary Extraction Balloons

- no allergic reactions (key word Latex allergy)
- extremely good visibility under fluoroscopy
- radiopaque markers with higher density value
- double- and triple lumen types
- huge variety on sizes and diameter
- single packed sterile



THE NEW GENERATION!

CE0483

STERILE EO

LATEX FREE



Order No.	Catheter Ø	Balloon Ø	length	Air Volume	Lumen	accepts guide wire
N135000	5 Fr.	9 mm	200 cm	1.5 ml	2	.025"
N135005	5 Fr.	12 mm	200 cm	2.8 ml	2	.025"
N135010	7 Fr.	14 mm	200 cm	2.8 ml	2	.035"
N135011	5 / 7 Fr.	14 mm	200 cm	2.5 ml	2	.035"
N135020	7 Fr.	18 mm	200 cm	5.0 ml	2	.035"
N135030	7 Fr.	16 mm	200 cm	5.0 ml	3	.035"
N135031	7 Fr.	14 mm	200 cm	4.0 ml	3	.035"
N135032	7 Fr.	12 mm	200 cm	2.5 ml	3	.035"
N135033	7 Fr.	9 mm	200 cm	1.8 ml	3	.035"

For more information about the products please contact your distributor.

09 / 2009



ENDO-FLEX®

ENDO-FLEX GmbH • Alte Hünxer Straße 115 • 46562 Voerde • Germany

☎ +49 (0)281 / 9 44 00 - 0
Fax. +49 (0)281 / 9 44 00 - 11
e-mail: export@endo-flex.de
<http://www.endo-flex.de>