



# ENDO-NEWS

## Extraction Balloon - 3 progressive stages -

- double or triple lumen type
- 7 Fr. sheath diameter
- fitting to working channel 2.8 mm
- for guide wire .035 inch
- improved radiopaque markers
- three different syringes included



**LATEX FREE!**



STERILE EO

CE0483

Order No.	lumen	catheter-Ø	catheter length	inflated balloon-Ø		
N135010PRO	double	7 Fr.	200 cm	9 mm	13 mm	16 mm
N135030PRO	triple					

For more information about the products please contact your distributor.

11 / 2010