



# ENDO-NEWS

## EasyCollect Polyp- & Foreign Body Retrieval Device



- oval thermo-formed plastic bag
- bag size 15 and 35 mm
- 1.8 and 2.3 mm Ø outer diameter
- latex-free
- rotatable
- single packed sterile
- for single-use only

Order No.	bag size	tube diameter	total length	packing unit
NE1518-B15	15 mm	1.8 mm	160 cm	10 pieces
NE1522-G35	35 mm	2.3 mm	230 cm	10 pieces

For more information about the products please contact your distributor.

11 / 2011