



# ENDO-NEWS

## Disposable Hot Biopsy Forceps

for current supported removal of small polyps



- oval fenestrated cups
- outer coating for highly effective insulation
- international HF connector
- ergonomical plastic handle
- single packed sterile



STERILE EO



CE0483

Order No.	Diameter	Length	Package Unit
NE6122-G	2.3 mm	230 cm	10 pcs

For more information about the products please contact your distributor.

12 / 2009



# ENDO-FLEX®

ENDO-FLEX GmbH • Alte Hünxer Straße 115 • 46562 Voerde • Germany

☎ +49 (0)281 / 9 44 00 - 0  
 Fax. +49 (0)281 / 9 44 00 - 11  
 e-mail: [export@endo-flex.de](mailto:export@endo-flex.de)  
<http://www.endo-flex.de>