



ENDO-NEWS

Disposable Cleaning Brush Set

Brush Set consist of:

- **Brush 1:** Gastroscope / Colonoscope cleaning brush
- **Brush 2:** Valve / Port brush

- double-sided cleaning brushes -
- perfectly rounded brush heads -



Order No.	Set components	Brush Dia.	Tube Dia.	Total Length	Packing Unit
NEBS41221-G-SET01	Brush 1	7 mm	1.8 mm	2.300 mm	25 pcs.
	Brush 2	5 and 12 mm	---	40 and 60 mm	

For more information about the products please contact your distributor.

06 / 2011