



Medi-Globe®

A MEDI-GLOBE COMPANY

Colonoscopy Accessories

Single Use



Polyp Trap™

For collection of polyps and specimens

Easy endoscope and suction device connectivity

Colonoscopy Shorts

Universal size

Colon Decompression Catheter

Radiopaque catheter with six sideports

Efficient decompression of toxic megacolon
and pseudo-obstructions

Colonoscopy Accessories – Single Use



Polyp Trap™

For polyp collection after excision of polyps and specimens for histological examination; sealable for transportation; including laboratory sticker sheet per box

Catalogue No.	Description
100132	Four clear, visible marked chambers and four open areas for straight suction. Pre-assembled highly flexible connection tubes for quick and easy installation.



Colonoscopy Shorts

With slit on backside; for protection of patient intimacy

Catalogue No.	Description
GAC-40-01-000	Bacteriostatic, water resistant and air permeable, universal size with slit on backside



Colon Decompression Catheter

Funnel connector; radiopaque tube

Catalogue No.	Description	Catheter Diameter (French)	Total Length (cm)	Compatible with Guide Wires
95016175	Radiopaque Tube with 6 Distal Sideports	16	175	.035" – .089"
95018175	Radiopaque Tube with 6 Distal Sideports	18	175	.035" – .089"



Medi-Globe GmbH

Medi-Globe-Straße 1–5 • D-83101 Rohrdorf OT Achenmühle
 Tel.: +49 (0) 8032 973-379 • Fax: +49 (0) 8032 973-392
 Email: sales@medi-globe.de • www.medi-globe.de