

# Stone Extraction Baskets

Single Use



L-GEB-SU-00-EN / LL-0032 Vers. 7.0 (12/2020)

## **Kink-resistant Nitinol Version**

Controlled radial deployment in the bile duct

## **45° Injection Port & Integrated Sealing**

User-friendly injection of contrast medium

## **Detachable Instrument Holder & Outer Sheath**

Applicable to emergency lithotripsy spiral and StoneBuster handle in emergency cases



### **Medi-Globe GmbH**

Medi-Globe-Straße 1–5 • D-83101 Rohrdorf OT Achenmühle  
Tel.: +49 (0) 8032 973-379 • Fax: +49 (0) 8032 973-392  
Email: [sales@medi-globe.de](mailto:sales@medi-globe.de) • [www.medi-globe.de](http://www.medi-globe.de)

# Stone Extraction Baskets – Single Use


# Stone Extraction Baskets – Single Use



## Stone Extraction Basket – 4 Braided Wires

Suitable for emergency lithotripsy system

Catalogue No.	Basket Shape	Tube Diameter (mm)	Distal Tip Diameter	Basket Width in Opened Condition (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GEB-11-19-420	4-Wire Basket	1.9	1.8	20	180	2.8
GEB-11-19-425	4-Wire Basket	1.9	1.8	25	180	2.8
GEB-11-19-430	4-Wire Basket	1.9	1.8	30	180	2.8
GEB-11-25-420	4-Wire Basket	2.35	2.6	20	180	2.8
GEB-11-25-425	4-Wire Basket	2.35	2.6	25	180	2.8
GEB-11-25-430	4-Wire Basket	2.35	2.6	30	180	2.8



## Stone Extraction Basket – 8 Monofilament Nitinol Wires

Not suitable for mechanical lithotripsy


Catalogue No.	Basket Shape	Tube Diameter (mm)	Distal Tip Diameter	Basket Width in Opened Condition (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GEB-13-19-820	8-Wire Dormia	1.9	1.8	20	180	2.8
GEB-13-19-830	8-Wire Dormia	1.9	1.8	30	180	2.8



## Stone Extraction Basket – 6 Braided Wires

Suitable for emergency lithotripsy system

Catalogue No.	Basket Shape	Tube Diameter (mm)	Distal Tip Diameter	Basket Width in Opened Condition (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GEB-11-19-620	6-Wire Basket	1.9	1.8	20	180	2.8
GEB-11-19-625	6-Wire Basket	1.9	1.8	25	180	2.8
GEB-11-19-630	6-Wire Basket	1.9	1.8	30	180	2.8
GEB-11-25-620	6-Wire Basket	2.35	2.5	20	180	2.8
GEB-11-25-625	6-Wire Basket	2.35	2.5	25	180	2.8
GEB-11-25-630	6-Wire Basket	2.35	2.5	30	180	2.8



## Stone Extraction Basket – 8 Monofilament Nitinol Wires

Compatible with guide wire .035", not suitable for mechanical lithotripsy

Catalogue No.	Basket Shape	Tube Diameter (mm)	Distal Tip Diameter	Basket Width in Opened Condition (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GEB-14-27-820	8-Wire Dormia	2.7	0.95	20	180	2.8
GEB-14-27-830	8-Wire Dormia	2.7	0.95	30	180	2.8