

RotaCut® Sphincterotomes

Single Use



L-GSP-SU-00-EN / LL-0067 Vers. 10.0 (01/2020)

Rotatable Catheter Tip

Precise positioning

2- and 3-lumen Versions

Use of contrast medium while guide wire is in place

Versions with Monofilament or Semi-Insulated Cutting Wire

Safe and controlled incision

Color Marking with Black Rings

Macroscopically controlled depth of
cutting wire penetration




Medi-Globe GmbH

Medi-Globe-Straße 1–5 • D-83101 Rohrdorf OT Achenmühle
Tel.: +49 (0) 8032 973-379 • Fax: +49 (0) 8032 973-392
Email: sales@medi-globe.de • www.medi-globe.de

RotaCut® Sphincterotomes – Single Use

RotaCut® 1-Lumen Sphincterotomes – Single Use

Sphincterotome Characteristics



RotaCut®: with multifilament cutting wire RotaCut®-M: with monofilament cutting wire Premier Cut: with multifilament cutting wire that is half proximally insulated (semi-insulated) S.O.S. Sphincterotomes: „Speed Option Set“, consisting of sphincterotome with integrated, preloaded .035“ „hydrophilic Nitinol guidewire in a length of 450 cm



RotaCut® – Sphincterotome with Braided Cutting Wire

Length 200 cm; colored ring markers

Catalogue-No.	Description	Tip Length (mm)	Cutting Wire Length (mm)	Outer Sheath Diameter (mm)	Min. Accessory Channel (mm)
GSP-01-17-020	Standard Tip	5	20	1.7	2.0



Standard Tip



Needle Knife Sphincterotome

Colored ring markers

Catalogue-No.	Needle Projection	Needle Stop at (mm)	Length (cm)	Outer Sheath Diameter (mm)	Min. Accessory Channel (mm)
GSP-11-17-005	Pointed Type	5	200	1.7	2.0
GSP-11-17-010	Pointed Type	10	200	1.7	2.0
GSP-11-17-015	Pointed Type	15	200	1.7	2.0



Pointed Type



RotaCut® 2-Lumen Sphincterotomes – Single Use



Tapered Tip

RotaCut® – with Braided Cutting Wire

Length 200 cm; compatible with guide wire .035" (0.89 mm); colored ring markers

Catalogue-No.	Description	Tip Length (mm)	Cutting Wire Length (mm)	Outer Sheath Diameter (mm)	Min. Accessory Channel (mm)
GSP-21-20-020	Tapered Tip	5	20	2.55/2.1	2.8
GSP-21-20-030	Tapered Tip	5	30	2.55/2.1	2.8
GSP-24-20-020	Tapered Tip	10 (long nose)	20	2.55/2.1	2.8
GSP-24-20-030	Tapered Tip	10 (long nose)	30	2.55/2.1	2.8
GSP-21-17-020	Tapered Tip	5	20	1.7	2.0
GSP-21-17-030	Tapered Tip	5	30	1.7	2.0



Billroth II-Type

RotaCut® – Billroth II-Type with Multifilament Cutting Wire

Length 200 cm; compatible with guide wire .035" (0.89 mm); colored ring markers

Catalogue-No.	Description	Tip Length (mm)	Cutting Wire Length (mm)	Outer Sheath Diameter (mm)	Min. Accessory Channel (mm)
GSP-34-20-020	Billroth II-Type	–	Variable	2.55/2.1	2.8



Tapered Tip

Premier Cut – with Braided and Semi-Insulated Cutting Wire

Length 200 cm; compatible with guide wire .035" (0.89 mm); colored ring markers

Catalogue-No.	Description	Tip Length (mm)	Cutting Wire Length (mm)	Outer Sheath Diameter (mm)	Min. Accessory Channel (mm)
GSP-25-20-020	Tapered Tip	5	20	2.55/2.1	2.8
GSP-25-20-030	Tapered Tip	5	30	2.55/2.1	2.8



RotaCut® 3-Lumen Sphincterotomes – Single Use




Tapered Tip

RotaCut® – with Braided Cutting Wire

Length 200 cm; compatible with guide wire .035" (0.89 mm); colored ring markers

Catalogue-No.	Description	Tip Length (mm)	Cutting Wire Length (mm)	Outer Sheath Diameter (mm)	Min. Accessory Channel (mm)
GSP-31-20-020	Tapered Tip	5	20	2.55/2.1	2.8
GSP-31-20-025	Tapered Tip	5	25	2.55/2.1	2.8
GSP-31-20-030	Tapered Tip	5	30	2.55/2.1	2.8



Tapered Tip

RotaCut®-M – with Monofilament Cutting Wire

Length 200 cm; compatible with guide wire .035" (0.89 mm); colored ring markers

Catalogue-No.	Description	Tip Length (mm)	Cutting Wire Length (mm)	Outer Sheath Diameter (mm)	Min. Accessory Channel (mm)
GSP-32-20-020	Tapered Tip	5	20	2.55/2.1	2.8
GSP-32-20-025	Tapered Tip	5	25	2.55/2.1	2.8
GSP-32-20-030	Tapered Tip	5	30	2.55/2.1	2.8



Tapered Tip

Premier Cut – with Braided and Semi-Insulated Cutting Wire

Length 200 cm; compatible with guide wire .035" (0.89 mm); colored ring markers

Catalogue-No.	Description	Tip Length (mm)	Cutting Wire Length (mm)	Outer Sheath Diameter (mm)	Min. Accessory Channel (mm)
GSP-27-20-030	Tapered Tip	5	30	2.55/2.1	2.8



RotaCut® 3-Lumen Sphincterotomes – Single Use



Pointed Type


Needle Knife Sphincterotome

Compatible with guide wire .035" (0.89mm); colored ring markers

Catalogue-No.	Needle Projection	Needle Stop at (mm)	Length (cm)	Outer Sheath Diameter (mm)	Min. Accessory Channel (mm)
GSP-15-25-005	Pointed Type	5	200	2.55/2.1	2.8
GSP-15-25-010	Pointed Type	10	200	2.55/2.1	2.8
GSP-15-25-015	Pointed Type	15	200	2.55/2.1	2.8



Pre-Loaded Sphincterotome Sets – Single Use



Tapered Tip

S.O.S. Sphincterotome Set – with Pre-Loaded Guide Wire

„Speed-Option Set“: 2-lumen sphincterotome with integrated pre-loaded and hydrophilically Teflon-coated .035" Nitinol-guide wire, length 450 cm (GGW-03-35-450); colored ring markers

Catalogue-No.	Description	Tip Length (mm)	Cutting Wire Length (mm)	Outer Sheath Diameter (mm)	Min. Accessory Channel (mm)
GSP-21-20-031	Tapered Tip Braided Cutting Wire	5	30	2.55/2.1	2.8
GSP-25-20-031	Tapered Tip Braided Cutting Wire Semi-Insulated	5	30	2.55/2.1	2.8



Tapered Tip

S.O.S Sphincterotome Set – with Pre-Loaded Guide Wire

Speed Option Set“: 3-lumen sphincterotome with integrated pre-loaded and hydrophilically Teflon-coated .035" Nitinol-guide wire, length 450 cm (GGW-03-35-450); colored ring markers

Catalogue-No.	Description	Tip Length (mm)	Cutting Wire Length (mm)	Outer Sheath Diameter (mm)	Min. Accessory Channel (mm)
GSP-32-20-031	Tapered Tip Monofilament Cutting Wire	5	30	2.55/2.1	2.8