

# Sphincterotomes

Single Use



LEF-GSP-SU-00-EN / LL-0336 Vers. 1.0 (07/2020)

## Spatial Tip Orientation

Precise positioning

## 2- and 3-Lumen Versions

Contrast medium injection while guide wire is in place

## Color Marking


Indicates macroscopically controlled cannulation depth



### ENDO-FLEX GmbH

Alte Hünxer Str. 115 • 46562 Voerde • Germany  
Tel.: +49 281 94400-0 • Fax: +49 281 94400-11  
Email: sales@endo-flex.de • www.endo-flex.de

## Sphincterotome Characteristics



## Sphincterotome without Guide Wire Access

Length 200 cm; colored ring markers



Pre-Cut-Type



Pre-Cut-Type with Tapered Tip



Tapered Tip

Catalogue No.	Description	Tip Length (mm)	Cutting Wire Length (mm)	Outer Sheath Diameter (mm)	Min. Accessory Channel (mm)
OE1031821	Pre-Cut-Type	-	20	1.8	2.2
OE1031831	Pre-Cut-Type	-	30	1.8	2.2
OE1051821	Tapered Tip, Strengthened	5	20	1.8	2.2
OE1051830	Tapered Tip, Strengthened	5	30	1.8	2.2
OE1051831	Tapered Tip, Strengthened	5	30	1.8	2.2
OE1051815PC	Pre-Cut-Type with Tapered Tip, Strengthened	-	15	1.8	2.2
OE1051821PC	Pre-Cut-Type with Tapered Tip, Strengthened	-	20	1.8	2.2
OE1051831PC	Pre-Cut-Type with Tapered Tip, Strengthened	-	30	1.8	2.2



## Sphincterotome compatible with Guide Wire .035"

Length 200 cm; colored ring markers



Pre-Cut-Type



Tapered Tip

Catalogue No.	Description	Tip Length (mm)	Cutting Wire Length (mm)	Outer Sheath Diameter (mm)	Min. Accessory Channel (mm)
OE1031830GW	Pre-Cut-Type	-	30	2.45/1.85	2.8
OE1051820GW	Tapered Tip	5	20	2.45/1.85	2.8
OE1051830GW	Tapered Tip	5	30	2.45/1.85	2.8
OE1051825GW11M	Tapered Tip, Monofilament Cutting Wire, Cutting Direction 11 o'clock	5	25	2.45/1.85	2.8



# 1-Lumen Sphincterotomes



B-II-Type

## B-II-Sphincterotome compatible with Guide Wire .035"

Catalogue No.	Description	Cutting Wire Length (mm)	Length (cm)	Outer Sheath Diameter (mm)	Min. Accessory Channel (mm)
OE11222FD	B-II-Type	15	200	2.8	3.0
OE11222FD-280*	B-II-Type	15	280	2.8	3.0



\* Packing unit 1 pc.



Pointed Type Blunt Type

## Needle Knife Sphincterotome without Guide Wire Access

Colored ring markers

Catalogue No.	Needle Projection	Needle Stopp at (mm)	Length (cm)	Outer Sheath Diameter (mm)	Min. Accessory Channel (mm)
OE11018N3	Pointed Type	3	200	1.8	2.0
OE11018N6	Pointed Type	6	200	1.8	2.0
OE11018N3-280	Pointed Type	3	280	1.8	2.0
OE11018SP3	Blunt Type	3	200	1.8	2.0
OE11018SP6	Blunt Type	6	200	1.8	2.0



# 2-Lumen Sphincterotomes



Standard Tip Pre-Cut-Type

## Sphincterotome compatible with Guide Wire .035"

Colored ring markers

Catalogue No.	Description	Tip Length (mm)	Cutting Wire Length (mm)	Length (cm)	Outer Sheath Diameter (mm)	Min. Accessory Channel (mm)
OE1002220DL	Standard Tip	5	20	200	2.55/2.1	2.8
OE1002225DL	Standard Tip	5	25	200	2.55/2.1	2.8
OE1002230DL	Standard Tip	5	30	200	2.55/2.1	2.8
OE1022220DL	Pre-Cut-Type	-	20	200	2.55/2.1	2.8
OE1022225DL	Pre-Cut-Type	-	25	200	2.55/2.1	2.8
OE1022230DL	Pre-Cut-Type	-	30	200	2.55/2.1	2.8



Smart Tip

## Sphincterotome "Smart Tip" compatible with Guide Wire .035"

Colored ring markers

Catalogue No.	Description	Tip Length (mm)	Cutting Wire Length (mm)	Length (cm)	Outer Sheath Diameter (mm)	Min. Accessory Channel (mm)
OE1042220DL	Tapered Tip	5	20	200	2.55/2.1	2.8
OE1042225DL	Tapered Tip	5	25	200	2.55/2.1	2.8
OE1042230DL	Tapered Tip	5	30	200	2.55/2.1	2.8
OE1042220DL-2	Tapered Tip	2	20	200	2.55/2.1	2.8
OE1042230DL-280*	Tapered Tip	5	30	280	2.55/2.1	2.8



\* Packing unit 1 pc.

# 2-Lumen Sphincterotomes



Smart Tip

## Preloaded Guide Wire Sphincterotome incl. HI-FLEX-Wire 21735450

Colored ring markers

Catalogue No.	Description	Tip Length (mm)	Cutting Wire Length (mm)	Length (cm)	Outer Sheath Diameter (mm)	Min. Accessory Channel (mm)
OE1042225DL-SET01	Smart Tip	5	25	200	2.55/2.1	2.8
OE1042230DL-SET01	Smart Tip	5	30	200	2.55/2.1	2.8



Smart Cut

## Sphincterotome "Smart Cut" compatible with Guide Wire .035"

Colored ring markers; partially insulated cutting wire

Catalogue No.	Description	Tip Length (mm)	Cutting Wire Length (mm)	Length (cm)	Outer Sheath Diameter (mm)	Min. Accessory Channel (mm)
OE1042220DLC	Smart Tip	5	20/10	200	2.55/2.1	2.8
OE1042225DLC	Smart Tip	5	25/12	200	2.55/2.1	2.8
OE1042230DLC	Smart Tip	5	30/15	200	2.55/2.1	2.8
OE1042230DLC1	Tapered Tip	5	30/15	200	2.55/2.1	2.8



# 2-Lumen Sphincterotomes



Pointed Type

## Needle Knife Sphincterotome compatible with Guide Wire .035"

Colored ring markers

Catalogue No.	Needle Projection	Needle Stopp at (mm)	Length (cm)	Outer Sheath Diameter (mm)	Min. Accessory Channel (mm)
OE11018N-DL	Pointed Type	6	200	2.55/2.1	2.8
OE11018N3-DL	Pointed Type	3	200	2.55/2.1	2.8



# 3-Lumen Sphincterotomes



Standard Tip



Pre-Cut-Type

## Sphincterotome compatible with Guide Wire .035"

Length 200 cm; colored ring markers

Catalogue No.	Description	Tip Length (mm)	Cutting Wire Length (mm)	Outer Sheath Diameter (mm)	Min. Accessory Channel (mm)
OE1002220TRL	Standard Tip	5	20	2.55/2.1	2.8
OE1002225TRL	Standard Tip	5	25	2.55/2.1	2.8
OE1002230TRL	Standard Tip	5	30	2.55/2.1	2.8
OE1022220TRL	Pre-Cut-Type	-	20	2.55/2.1	2.8
OE1022225TRL	Pre-Cut-Type	-	25	2.55/2.1	2.8
OE1022230TRL	Pre-Cut-Type	-	30	2.55/2.1	2.8
OE1002220TRL-M	Standard Tip, Monofilament Cutting Wire	5	20	2.55/2.1	2.8
OE1002225TRL-M	Standard Tip, Monofilament Cutting Wire	5	25	2.55/2.1	2.8
OE1002230TRL-M	Standard Tip, Monofilament Cutting Wire	5	30	2.55/2.1	2.8



# 3-Lumen Sphincterotomes



Tapered Tip

## Sphincterotome "Tapered Tip" compatible with Guide Wire .035"

Length 200 cm; colored ring markers

Catalogue No.	Description	Tip Length (mm)	Cutting Wire Length (mm)	Outer Sheath Diameter (mm)	Min. Accessory Channel (mm)
OE1042220TRL	Tapered Tip	5	20	2.55/2.1	2.8
OE1042225TRL	Tapered Tip	5	25	2.55/2.1	2.8
OE1042230TRL	Tapered Tip	5	30	2.55/2.1	2.8
OE1042220TRL-3	Tapered Tip	3	20	2.55/2.1	2.8
OE1042220TRL-M	Tapered Tip, Monofilament Cutting Wire	5	20	2.55/2.1	2.8
OE1042225TRL-M	Tapered Tip, Monofilament Cutting Wire	5	25	2.55/2.1	2.8
OE1042230TRL-M	Tapered Tip, Monofilament Cutting Wire	5	30	2.55/2.1	2.8



# 3-Lumen Sphincterotomes



Smart Tip

## Preloaded Guide Wire Sphincterotome incl. HI-FLEX-Wire 21735450

Length 200cm; colored ring markers

Catalogue No.	Description	Tip Length (mm)	Cutting Wire Length (mm)	Outer Sheath Diameter (mm)	Min. Accessory Channel (mm)
OE1042220TRL-SET01	Tapered Tip	5	20	2.55/2.1	2.8
OE1042225TRL-SET01	Tapered Tip	5	25	2.55/2.1	2.8
OE1042225TRL-M-SET01	Tapered Tip , Monofilament Cutting Wire	5	25	2.55/2.1	2.8
OE1042230TRL-M-SET01	Tapered Tip, Monofilament Cutting Wire	5	30	2.55/2.1	2.8



Smart Cut

## Sphincterotome "Smart Cut" compatible with Guide Wire .035"

Length 200cm; colored ring markers; partially insulated cutting wire

Catalogue No.	Description	Tip Length (mm)	Cutting Wire Length (mm)	Outer Sheath Diameter (mm)	Min. Accessory Channel (mm)
OE1042220TRLC	Tapered Tip	5	20/10	2.55/2.1	2.8
OE1042225TRLC	Tapered Tip	5	25/12	2.55/2.1	2.8
OE1042230TRLC	Tapered Tip	5	30/15	2.55/2.1	2.8



# 3-Lumen Sphincterotomes



Pointed Type

## Needle Knife Sphincterotome compatible with Guide Wire .035"

Colored ring markers

Catalogue No.	Needle Projection	Needle Stopp at (mm)	Length (cm)	Outer Sheath Diameter (mm)	Min. Accessory Channel (mm)
OE11018N-TRL	Pointed Type	6	200	2.55/2.1	2.8
OE11018N4-TRL	Pointed Type	4	200	2.55/2.1	2.8

