

# Cysto-Gastro Sets

Single Use



## **For Ultrasound-Guided Endoscopies**

For transgastric and transduodenal puncture,  
and for pancreatic pseudocysts

## **6 French Version incl. Cystotome**

Designed for use in small working channels

## **„2 in 1“ – 8.5 and 10 French Version incl. Cystotome & Puncture Needle**

For placement of drainage stents

# Cysto-Gastro-Sets – Single Use



## Cysto-Gastro-Set<sup>\*,\*\*</sup>

Includes needle knife; outer catheter with diathermic tip

Catalogue No.	Description	Working Length (cm)	Catheter Outer Diameter (French)	Compatible Guide Wire (inch)
E30008341	2 HF-Connectors	180	8.5	.025
E30010340	2 HF-Connectors	180	10	.035



## Cysto-Gastro-Set/Cystotome<sup>\*\*,\*\*\*</sup>

With diathermic tip

Catalogue No.	Description	Working Length (cm)	Catheter Outer Diameter (French)	Compatible Guide Wire (inch)
E30006345	1 HF-Connector	180	6	.035

### \* Fr. Type 10 and 8.5

After the initial access with the implemented needle knife, the incision can be enlarged with the diathermic metal tip on the Fr. 10 or Fr. 8.5 outer sheath. The sets are equipped with two HF connectors, one each for the needle knife and the diathermic tip on the outer catheter.

### \*\* Drainage

After access to the identified pseudocyst appropriate stents can be placed over the guide wire to facilitate drainage from the cyst to either the stomach or the duodenum.

### \*\*\* Fr. Type 6

Specifically designed for smaller working channels, the Fr. 6 Cystotome has a diathermic metal tip on the Fr. 6 catheter and one HF connector only. It is advanced over a .035" guide wire to the transgastric or transduodenal wall directly.

LEF-CGS-SU-00-EN / LL-0354 Vers. 2.0 (08/2020)



### ENDO-FLEX GmbH

Alte Hünxer Str. 115 • 46562 Voerde • Germany  
 Tel.: +49 281 94400-0 • Fax: +49 281 94400-11  
 Email: sales@endo-flex.de • www.endo-flex.de