

# Polypectomy / Mucosectomy Snares

Single Use




LEF-GPS-SU-00-EN / LL-0363 Vers. 1.0 (08/2020)

**Versions with 360° Rotation**  
Easier polyp snaring

**Monofilament Loop Versions**  
Precise and smooth cutting

## Rotatable Polypectomy Snares – Single Use



Hexagonal Loop



### Polypectomy Snare

Hexagonal loop, braided wire, rotatable

Catalogue No.	Loop Shape	Loop Diameter (mm)	Sheath Diameter (mm)	Min. Accessory Channel (mm)	Working Length (cm)
NOE372214-G	Hexagonal	15	2.35	2.8	230
NOE372216-G	Hexagonal	25	2.35	2.8	230
NOE372217-G	Hexagonal	35	2.35	2.8	230



Asymmetrical Loop



### Polypectomy Snare

Asymmetrical loop, braided wire, rotatable

Catalogue No.	Loop Shape	Loop Diameter (mm)	Sheath Diameter (mm)	Min. Accessory Channel (mm)	Working Length (cm)
NOE3222-G	Asymmetrical	15	2.35	2.8	230
NOE3322-G	Asymmetrical	25	2.35	2.8	230



**ENDO-FLEX GmbH**

Alte Hünxer Str. 115 • 46562 Voerde • Germany  
 Tel.: +49 281 94400-0 • Fax: +49 281 94400-11  
 Email: sales@endo-flex.de • www.endo-flex.de

# Rotatable Polypectomy Snares – Single Use



Oval Loop

## Polypectomy Snare

Oval loop, braided wire, rotatable

Catalogue No.	Loop Shape	Loop Diameter (mm)	Sheath Diameter (mm)	Min. Accessory Channel (mm)	Working Length (cm)
NOE341816-C	Oval	25	1.8	2.0	180
NOE342213-G	Oval	10	2.35	2.8	230
NOE342214-C	Oval	15	2.35	2.8	180
NOE342214-G	Oval	15	2.35	2.8	230
NOE342214-M	Oval	15	2.35	2.8	280
NOE342216-C	Oval	25	2.35	2.8	180
NOE342216-G	Oval	25	2.35	2.8	230
NOE342217-C	Oval	35	2.35	2.8	180
NOE342217-G	Oval	35	2.35	2.8	230
NOE342217-M	Oval	35	2.35	2.8	280



Oval Loop

## Polypectomy Snare

Oval loop, monofilament wire, rotatable

Catalogue No.	Loop Shape	Loop Diameter (mm)	Sheath Diameter (mm)	Min. Accessory Channel (mm)	Working Length (cm)
NOE342212M-G	Oval	6	2.35	2.8	230
NOE342213M-G	Oval	10	2.35	2.8	230
NOE342214M-C	Oval	15	2.35	2.8	180
NOE342214M-G	Oval	15	2.35	2.8	230
NOE342216M-C	Oval	25	2.35	2.8	180
NOE342216M-G	Oval	25	2.35	2.8	230
NOE342217M-G	Oval	35	2.35	2.8	230



# Polypectomy Snares – Single Use



Oval Loop

## Polypectomy Snare

Oval loop, braided wire

Catalogue No.	Loop Shape	Loop Diameter (mm)	Sheath Diameter (mm)	Min. Accessory Channel (mm)	Working Length (cm)
OE342215-G	Oval	20	2.3	2.8	230
OE342216-G	Oval	25	2.3	2.8	230
OE342217-C	Oval	35	2.3	2.8	180
OE342217-G	Oval	35	2.3	2.8	230



Oval Loop

## Polypectomy Snare

Oval loop, monofilament wire

Catalogue No.	Loop Shape	Loop Diameter (mm)	Sheath Diameter (mm)	Min. Accessory Channel (mm)	Working Length (cm)
OE342214M-G	Oval	15	2.3	2.8	230



Asymmetrical Loop

## Polypectomy Snare

Asymmetrical loop, braided wire

Catalogue No.	Loop Shape	Loop Diameter (mm)	Sheath Diameter (mm)	Min. Accessory Channel (mm)	Working Length (cm)
OE3322-G	Asymmetrical	25	2.3	2.8	230



# Mucosectomy Snares – Single Use



Oval Loop

## Mucosectomy Snare

Oval loop, braided wire

Catalogue No.	Loop Shape	Loop Diameter (mm)	Sheath Diameter (mm)	Min. Accessory Channel (mm)	Working Length (cm)
OE342246-G	Oval	25	2.3	2.8	230

